Progress of laser cooling of ¹²C³⁺ ions at the CSRe

H. B. Wang *^{, † 1}, X. Ma * ², W. Q. Wen*, Z. K. Huang^{*, †}, D. C. Zhang^{*}, B. Hai^{*, †}, X. L. Zhu^{*}, D. M. Zhao*, J. Li*, X. M. Ma*, T. L. Yan*, R. S. Mao*, T. C. Zhao*, J. X. Wu*, J. C. Yang*, Y. J. Yuan*, J. W. Xia*, M. Loeser[&], M. Siebold[&], U. Schramm[&], O. Boine-Frankenheim[§], L. Eidam[§], D. Winters[§], G. Birkl[#], B. Rein[#], Th. Walther[#], and M. Bussmann[&]

* Institute of Modern Physics, Chinese Academy of Sciences, 730000 Lanzhou, China University of Chinese Academy of Sciences, 100049 Beijing, China [§]GSI Helmholtzzentrum für Schwerionenforschung GmbH, 64291 Darmstadt, Germany [#] Institut für Angewandte Physik, Technische Universität Darmstadt, 64289 Darmstadt, Germany

& Helmholtz-Zentrum Dresden-Rossendorf, 01328 Dresden, Germany

Synopsis We have performed a test run for laser cooling experiments with ${}^{12}C^{3+}$ ion beams at an energy of 122 MeV/u at the CSRe with a pulsed laser. During this beamtime a lot of progress have been made. This was the first time we could successfully separate the ¹²C³⁺ ions and ¹⁶O⁴⁺ ions in the Schottky spectrum with the help of electron cooling. And our newly installed CPM detector worked well during the experiment. We tried to see effects from the interaction of the pulsed laser light with the stored ${}^{12}C^{3+}$ ion beams, but no cooling effects could yet be observed.


Laser cooling is the most promising method to achieve high phase-space densities or even crystalline beams for relativistic heavy ion at storage rings [1]. In order to achieve laser cooling at the CSRe in IMP, we performed a test run for a laser cooling experiment with ${}^{12}C^{3+}$ ion beams at an energy of 122 MeV/u at the CSRe with a pulsed laser in September 2014.

During this beamtime we have made a lot of progress and managed to get a large part of the experiment under control. Firstly, we could successfully separate ¹²C³⁺ ions and ¹⁶O⁴⁺ ions in the Schottky spectrum with the help of electron cooling for the first time (Figure 1). Secondly, the new CPM (Channeltron photomultiplier) detector which was used to detect the fluorescence from the (laser excited) $2p_{1/2}$ and the $2p_{3/2}$ states was tested systematically. With the laser off, the CPM detector was also used to determine the lifetime by recording fluorescence due to charge-changing collisions between ¹²C³⁺ ions and the residual gas. The lifetime was about 20 seconds, which is much longer than it was during the last beam time. Thirdly, we tried to see effects from laser cooling with the pulsed laser, and the beam and laser widths were measured using the scrapers and a wire scanner to make sure they were overlapped. Moreover, the dynamics of the electron-cooled and RFbunched ion beams were investigated systematically by longitudinal modulation of sinusoidal, square and barrier bucket waveforms.

¹ E-mail: <u>wanghanbing@impcas.ac.cn</u>

² E-mail: x.ma@impcas.ac.cn

⁰s 50s 12 (3+ 100s 150s With e-cooling 200s -50k -100k 0k 50k 100k f - 244.0104 M [Hz]



We tried to see effects from the interaction of the pulsed laser with the stored ${\rm ^{12}C^{3+}}$ ion beams. But unfortunately, we did not find any evidence for effects from the interaction in the Schottky spectrum or the count rate of the CPM detector. We attribute this to the background and the low repetition rate of the laser. Therefore, we will repeat the experiment in 2015 using a laser system with a higher repetition rate.

This work is supported by the NSFC No. 10921504 and the External Cooperation Program of the CAS, Grant No. GJHZ1305. We acknowledge Tobias Beck for his help for preparing the beam time.

References

[1] U. Schramm et al 2004 Progress in Particle and Nuclear Physics 53 583-67

Content from this work may be used under the terms of the Creative Commons Attribution 3.0 licence. Any further distribution (\mathbf{i}) (cc) of this work must maintain attribution to the author(s) and the title of the work, journal citation and DOI. Published under licence by IOP Publishing Ltd